

Advantages

- ✓ Direct, simple, and economic point-of-care determination of cholinesterase activities
- ✓ Whole blood testing
 - one drop of blood (10 µl) per test
- ✓ The result is available after only 5 minutes
- ✓ Comfortable handling via touchscreen
- ✓ Portable due to the use of Li-Ion battery technology
- ✓ Integrated barcode reader for patient IDs
- ✓ Export of data via USB

Presentation of the results



Manufacturer:
Dr. Franz Köhler Chemie GmbH
Werner-von-Siemens-Str. 14-28
64625 Bensheim, Germany
Tel. +49 6251 1083-0, Fax +49 6251 1083-146
info@koehler-chemie.de, www.koehler-chemie.de

DR. F. KÖHLER CHEMIE

Pharmaceuticals



LISA-CHE

Point-of-care measurement of cholinesterase activity



Portable point-of-care device for the activity determination of acetyl- and butyrylcholinesterase

DR. F. KÖHLER CHEMIE

Pharmaceuticals



Measuring ChE-activities – Why?

Since 1959, the measurement of cholinesterase activity is used in patients with organophosphate poisoning.

In recent years, it has been shown that the neurotransmitter acetylcholine, acting as a parasympathetic messenger, regulates the body's immune response via the cholinergic anti-inflammatory reflex. This reflex controls the release of proinflammatory cytokines and, thus, neuroinflammation.

Furthermore, it was shown that the activities of cholinesterases can serve as biomarkers, for example

- for the detection of a postoperative delirium
- in systemic inflammation
- as a predictor of the mortality in patients with extracorporeal membrane oxygenation

Measurement of cholinesterase activity for neuroinflammation and delirium

When to measure?

1. Pre-anesthetic interview (→ base value)
 2. Pre-operatively
 3. After admission to ICU
- Optional: perioperatively

Delirium detection^{1,2}

1. Hyperactive delirium

- a) Strong hyperactive type:
 - relatively easy to identify
- b) Moderate type:
 - Measure BChE activity
 - Methods for delirium detection (e.g. CAM-ICU)

2. Hypoactive delirium

- a) Moderate type:
 - Measure BChE activity
 - Methods for delirium detection (e.g. CAM-ICU)
- b) Strong hypoactive type:
 - Measure BChE activity


Interpretation of BChE values³

- BChE activity > 2000 U/l: Normal range
- BChE activity < 2000 U/l: Increased attention necessary
- BChE activity < 1500 U/l: Suspected neuroinflammation/delirium
 - Additionally, if possible:
 - Methods for delirium detection (e.g. CAM-ICU)

Source references:

- [1] Maldonado JR, Acute Brain Failure: Pathophysiology, Diagnosis, Management, and Sequelae of Delirium, Crit Care Clin. 2017 Jul;33(3):461-519.
- [2] Cerejeira J et al., Low preoperative plasma cholinesterase activity as a risk marker of post-operative delirium in elderly patients, Age Ageing. 2011 40(5):621-6.
- [3] Performance assessment by Dr. Franz Koehler Chemie GmbH

An example: AChE and BChE in a patient with hypoactive delirium



Barth, E. (2014). Anticholinium® (physostigmine) as a therapeutic option in ICU patients with postoperative delirium and vigilance deficiency with cognitive dysfunction. Poster from DIVI, Hamburg

The graph shows the AChE- and BChE activities of a patient with a hypoactive delirium, who was treated with physostigmine (**ANTICHOLIUM**®) because of altered ChE-activities. Within one hour after the treatment, no delirium could be observed. 24 hours after the treatment, the AChE- and BChE activities were in the normal range.